

Scalable instrument management system improves accuracy and eliminates errors

Unfortunately, inaccurate paper count sheets and tray errors frequently occur in surgical facilities because of illegible mark-ups from multiple users, product confusion and other human errors.

To help alleviate these costly mistakes and eliminate your vast inventory of paper count sheets, SourceMedical offers the SourcePlus Instrument Manager, powered by CENSIS®. Instrument Manager is an effective, easy to implement, affordable surgical instrument management solution that configures instrument sets to a precise standard and improves the accuracy of your tray content.

With SourcePlus Instrument Manager, your facility will experience:

- **Improved count sheet accessibility and tray accuracy – saving labor time**

The system is comprised of a database of information that allows technicians to simply scan a tray barcode and print the correct count sheet for accurate assembly and precise count sheet management.

Assembly and staff training are also made easier using the database of information and photographs of more than 40,000 commonly used surgical instruments, along with images of your customized tray configuration to assist with any discrepancies.

With a single electronic repository of information, the staff now has one central location for managing tray count sheets. Staff will no longer have to search through filing cabinets for the correct count sheet or make photocopies each time a revision occurs.

- **User-friendly technology working behind the scenes**

The SourcePlus Instrument Manager's touch-screen, menu-driven user-interface allows any staff member to perform the required tasks with minimal training. And, because the system is web-based, updates of the latest software versions will occur without costly downtime. As additional protection, your data is backed up regularly by an off-site server, allowing you and your staff to focus more on patient care and less on the technical details.

SourcePlus Products and Services: Solving problems and boosting performance

The SourcePlus suite includes a wide array of add-on modules, products and services that add value to our core information management systems; further enhancing facility management, clinical efficiency and profitability. With our integrated SourcePlus Products, your surgical facility can leverage the best tools available to ASCs and surgical hospitals.

Benefits include:

- Improve service to patients and physicians
- Optimize reimbursement and cash flow
- Improve clinical documentation
- Boost staff productivity
- Reduce staff and supply costs
- Improve decision-making
- Improve customer service

The SourcePlus Instrument Manager offers seven easy steps for your facility staff to simplify tray processing:

1. Receive tray as it exits the washer
2. Scan the tray barcode to print the count sheet
3. Assemble instruments in tray
4. Verify tray contents are accurate and complete
5. Build a sterilization load with assembled containers, wraps, etc.
6. Assign the load to a sterilizer via the bar code scanner
7. Perform interface with sterilizers to automatically capture sterilization records (optional)

Scalable Instrument Management Solution

SourcePlus Instrument Manager is the perfect solution for many surgical facilities. For larger and more complex organizations requiring more precise management at either the tray or individual instrument levels, SourcePlus Instrument Manager can be upgraded to more comprehensive Tray or Instrument Management solutions powered by W. These scalable management solutions provide increased tray accuracy and precise logistical inventory management of trays.



About SourceMedical

SourceMedical provides outpatient information solutions and services for ambulatory surgery centers, surgical hospitals and rehabilitation clinics nationwide. With a nearly 30-year track record and more than 4,100 satisfied customers, SourceMedical is the trusted source for innovative applications, in-depth industry expertise and unsurpassed customer service. The company's unique, end-to-end systems improve operational efficiency and cash flow while enabling healthcare facilities to deliver a higher standard of patient care.

To find out how SourceMedical can help you achieve your clinical and financial goals, call 800-719-1904 or visit www.sourcemed.net.

